

Lumify Tablet Compatibility



Tablet or Phone Name	Screen Size	Weight	Approximate Continuous Scan Time	Connectivity Options	Internal Storage	USB Type
Samsung Galaxy S7 Edge	5.5 in	157 g	2 h 15 min	Wi-Fi/Cellular	32 GB	B
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	260 g	2 h 30 min	Wi-Fi	32 GB	B
Samsung Galaxy Tab S2 9.7" (Only Model SM-T813)	9.7 in	395 g	3 h	Wi-Fi/LTE	32 GB	B
Samsung Galaxy View	18.4 in	3900 g	5 h	Wi-Fi/LTE	32 GB	B
Huawei Mediapad M3 (Model BTW-W09)	8.4 in	310 g	3 h 30 min	Wi-Fi	32 GB	B
Lenovo Yoga Tab 3 (Model YT3-850F)	8 in	834 g	4 h	Wi-Fi	16 GB	B
Lenovo Yoga Tab 3 (Model YT3-X50F)	10 in	1150 g	4 h	Wi-Fi	16 GB	B

Compatible with L12-4/C5-2 only (not S4-1)

Samsung Galaxy Tab S2 8" (Model SM-T710)	8 in	260 g	2 h 30 min	Wi-Fi	32 GB	B
Samsung Galaxy Tab S2 9.7" (Model SM-T810)	9.7 in	395 g	3 h	Wi-Fi/LTE	32 GB	B

Compatible, but discontinued by manufacturer

Google Nexus 6 (not Nexus 6P)	6 in	185 g	1 h 45 min	Wi-Fi/Cellular	32/64 GB	B
NVIDIA SHIELD Tablet (V1)	8 in	390 g	1 h 30 min	Wi-Fi	16 GB	B
Acer Iconia Tab 10 A3-A30	10 in	540 g	2 h 30 min	Wi-Fi	16 GB	B
Samsung Galaxy Tab Active	8 in	370 g	3 h 30 min	Wi-Fi	16 GB	B
Samsung Galaxy Tab S	8.5 in	204 g	2 h 15 min	Wi-Fi	16 GB	B
Google Nexus 9	9 in	425 g	2 h	Wi-Fi/LTE	16/32 GB	B
Nvidia Shield Tablet K1	8 in	350 g	2 h 30 min	Wi-Fi	16 GB	B
Samsung Galaxy Note5	5.7 in	171 g	3 h 15 min	Wi-Fi/Cellular	32/64/128 GB	B

Tablet compatibility is subject to change, and Philips reserves the right to make changes without notice.



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.

239034 Rev L July 2017